

ASTM LEVEL 3 MASK

PROCEDURE EAR LOOPS | BLUE

REF PG4-1273

ABOUT THE MASK

Based on ASTM F2100 standards, this mask provides the highest level of protection during medical procedures. Advanced fluid-resistant and filtration capabilities of the mask provide all-in-one protection for wound care, line changes and dental work. This mask is recommended for procedures with a high risk of fluid splash, when handling trauma or caring for ER and ICU patients.

PROTECTION STANDARDS



- BFE (%): ≥ 99
- PFE (%): ≥ 99
- ΔP (mmH₂O/cm²): < 5
- Fluid Resistance: 160mm Hg
- Flame Spread: Class 1

Variation in lot results exist. All results meet or exceed ASTM F2100 requirements.

\$17.95/box of 50

**BUY 6
FREE
GET 4**

Net \$10.77/box of 50

After free goods
Offer valid until Dec 31, 2022.

WOW

FEATURES & BENEFITS

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin.

WIDE NOSEPIECE

Constructed with high-quality aluminum to keep the mask in place across the nose and cheekbones.

SOFT MASK AND EAR LOOP MATERIAL

Carefully engineered and clinically preferred for extended periods of use without discomfort.



QUANTITIES

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500



P. (902) 454.8393 | 888.839.8393
F. (902) 454.8433 | 877.833.4905

Did you know you can get AIR MILES® REWARD MILES™ with every purchase from CAN-med Healthcare?

Contact us to find out how!



18005657553

roseanna.murphy@canmedhealthcare.com



©™ Trademarks of AM Royalties Limited Partnership used under license by LoyaltyOne, Co. and CAN-Med Healthcare